TRONG

CA

ONTRIN

Velopeds & Walkers Made for Walking

Made for Walking

We believe in quality, good engineering and well thought out design. We are constantly focusing on innovating and making products that you can perform and be active with, without compromising the actual need for support.

We offer you products that promote health, mental well-being, enjoyment of the outdoors and opportunities for new experiences and achievements.

We want you to enjoy the best product experience in the market of our field; the ultimate support in urban and outdoor areas.

We hold all our customers dear and we embrace your collaboration, feedback and thoughts as they make our vision stronger; to make the best walking aids in the world.

That is why we don't settle for anything less than excellence in everything we do. We are honest, hard working and have the courage to change what we can, where we can.

Welcome!

Stefan & Johan

Cover Photo:

Domenico and Alice from Italy enjoying the scenery of the beach.

SUBJECT TO NORMAL

User Photos:

All user photos featured in this brochure are photos of actual users, that we have received from our enthusiastic customers.

SWEDISH

DESIGN

Keen photographer Ann-Marie captures the beautiful spring view of Klingerfjärden in Medelpad, Sweden.

TRIONIC

Lisbeth on tember wal skär, Swede

Be active - stay active

Active lifestyle

If you lead an active lifestyle – a Trionic Veloped or Walker is your premier choice of walking aid. They are the state-of-the-art walking support. By choosing a Trionic you'll enjoy greater possibilities to be physically active outdoors, to stay healthy and fit, and to lead a truly active lifestyle.

Positive health effects

Daily physical activity is the single most important measure to stay healthy and live longer. By walking for a total of 30 minutes a day you will meet doctors' recommendations for physical activity.

We are happy to offer you the tools so that you can be active and stay active.

User reviews

Through the independent online services provided by Trustpilot, customers and users all over the world have rated Trionic's products as "Excellent", with an overall rating of 4.9 out of 5.0.

My Veloped has changed my life

"Since the operation to remove a tumour from my spinal cord, I can only walk with two crutches and this is a slow and painful process. I was finally told about rollators, but thought they looked feeble and unattractive. My wish was to walk on the beach, or on cobbles in town. After searching the internet I came across the Trionic Veloped, and with little hesitation ordered one. It is the best purchase I have ever made. I can now walk further and a little faster, with confidence and without the fear of falling. It has opened up my world again. I can't recommend the Trionc Veloped highly enough."

Cathryn Shilling, United Kingdom

WOW!!!!

"This is an amazing machine for walking. I've been a sporty person all my life and now need some help to live the life I want to live. Because they've put themselves in the place of people who will use it, everything works easily and intuitively – no need for complicated instructions. The build quality is excellent. Intelligent engineering by people who have thought it through in every detail. The bonus? – It's cool!"

Robert Heap, United Kingdom

Outstanding

"Although the Veloped Tour is an expensive investment, it is well worth it. It has given me mobility that no other rollator provides. I can now navigate uneven terrain, cobble stone streets, curbs and even pot holes with ease. I intend to take it with me to Europe for use on a river cruise from Budapest to Prague. Hopefully, it will enable me to keep up with the group during walking tours. I could not be happier with this amazing product.

Loeb Granoff, USA


Excellent Product, Excellent Service

"We do not have any good outdoor walkers down under in Australia. I found the Trionic Veloped on the internet and e-mailed the company. I received outstanding support and service from the founder and inventor Stefan. He arranged delivery to Australia for a very reasonable fee. The Veloped itself is engineering Gold! My client with a brain injury hasn't walked around his property in 2 years, until he got his Veloped. He is stoked and so am I! The design allows a lot of stability as the person almost walks "in" the frame as opposed to pushing it in front of him/her. This gives it a lot more sideways suppport. All the gadgets and storage options allow my client to participate in meaningful life roles. It is the Iphone of outdoor walkers."

Kyle Maloney,

Community rehabilitation physiotherapist, Australia





Features and Functions on all Velope



The Wheel Size Matters

The performance of a rollator or walker relies heavily on the size of the wheels. The bigger the wheel, the easier it is to pass obstacles and uneven ground, and the comfort level also increases with the wheel size. A bigger wheel simply performs better than a small wheel, especially outdoors.

Trionic offers three different wheel sizes: 9" (23 cm), 12" (31 cm) and 14" (36 cm).



Air-Pressured Tires

The Walker 12er/14er and the Velopeds are all equipped with air-pressured tires. There are several advantages with air-filled tires, and they simply offer you a softer ride and unbeatable comfort. As they do not cause vibration like solid tires do, you are less prone to suffer pain and discomfort in your hands, elbows and shoulders when you walk on uneven ground.

The Velopeds and the Walker 12er/14er have high-quality lightweight tires with puncture protection and reflective logos.



All-Terrain Brake

Thanks to the All-Terrain Brake system you can always rely on an optimal brake performance, regardless if you are walking in sunshine, pouring rain or deep snow. Instead of a conventional rollator brake, working directly on the tire surface, all Velopeds and Walkers are equipped with powerful and exclusive hub brakes that are fully protected from dirt, water and snow.

ds and Walkers



Trionic Brake Lever

Trionic has developed a brand new and user-friendly brake lever. The all-aluminium brake lever offers you a solid and high quality feel. When activating the parking brake the parking trigger is easy to operate, and thanks to its design it offers multiple ways to activate it.

The aluminium lever has a Kraton rubber surface that prevents your fingers from getting cold from touching metal when you are walking in winter conditions.



Detachable Wheels

All Velopeds & Walkers have detachable wheels. By using the quick-releases on the wheels, you can easily detach them in order to minimize the transport weight and dimensions of your Veloped or Walker.

As you can quickly change the wheels, this also allows you to use two different wheel sets; one set with proper off-road tires and another set with smooth city tires. All push-button quick-release axles are made of 100% stainless steel.



Folding Function

Trionic Veloped: Through a quick push of the button on the telescopic folding stay, you can fold your Veloped to pass narrow passages. To fully fold it for transport or storage, press the second button on the telescopic stay and fold the frame all the way.

Trionic Walker: The Trionic Walker has an X-shaped folding function, and by pulling the handle positioned on the seat you fold it in a second. The frame folds very compact, and transport and storage is a breeze.



Patent Pending:

Trionic Syncro Steering Avoid Shimmying Wheels

Shimmying and skewed swivel wheels are frequently mentioned by rollator users as being a definite problem area. When walking on uneven ground with a 4-wheeled rollator, sooner or later one of the swivel wheels loses contact with the ground, starts shimmying and regains its contact with the ground in a skewed position. As a result the rollator often comes to an abrupt stop, and there is a considerable risk of the user tripping or even falling over.

Trionic Syncro

The new Trionic Walker rollator is equipped with a synchronized steering feature (patent pending) that eliminates the problems with shimmying and skewed swivel wheels.

By connecting the swivel wheels with each other through a pulley/ timing belt arrangement, the two wheels steer in unison. If one swivel wheel suddenly loses contact with the ground, the other wheel will still keep it positioned in your desired walking direction.

An additional central pulley guarantees that you will enjoy the same type of steering that your car offers, called the "Ackermann Steering Geometry". When you turn, the inner wheel is allowed to turn at a slightly sharper angle than the outer wheel.

Thanks to its clever design, the function of the Trionic Syncro Steering doesn't get affected when you fold and unfold the rollator.

Patent:

Trionic Climbing Wheel

Overcoming Obstacles

Overcoming obstacles like curbs, roots and rocks is a very common problem among rollator users. If the obstacle is higher than 1-2 inches (depending on the wheel size), the rollator needs to be lifted in order to pass the obstacle. If the user is forced to lift the rollator, or to perform some other maneuver to overcome the obstacle, the rollator has completely lost its main function of offering support.

Trionic Climbing Wheel

The Veloped is equipped with the Trionic Climbing Wheel (Patent) that greatly increases the ability to overcome different kinds of obstacles when walking outdoors. By spreading the load over two wheels, slightly displaced in the traveling direction, the climbing sequence is divided into two steps.

The climbing capacity is enhanced even more through an angled link arm that uses the pushing force provided by the user. This gives the Veloped a unique feature, the capacity to climb a 5 inch high obstacle, for instance a high curb or rock. When overcoming a lower obstacle, for example a lowered curbstone at a zebra crossing, you only need a fraction of the pushing force a single wheel requires. The new Veloped with 14-inch wheels will climb such an obstacle without you noticing that it was even there.

The pivoting link arm also provides suspension. On uneven ground it constantly pivots and suspends the ride. The suspension turns an off-road walk with a Veloped into pure pleasure.

The Trionic Climbing Wheel of the Veloped 12er in action, climbing a 5 inch high curb.



Meller

The New Premium

9er • 12er • 14er

between different seat and basket colour styles. Learn more about the possible tire and grip options on pages 24-25.

Easy Folding

By pulling the seat handle upwards you quickly and easily fold the Walker. When folded the Walker stands upright by itself, and thanks to the X-folding mechanism it is very compact. To unfold the Walker, you pull the handles outwards, and lock the frame in walking mode by pressing the seat downwards.

All-Aluminium Design

The Walker is a true aluminium rollator, and it has no major frame, handle bar or folding cross parts made from plastic. The frame is made from a machine-formed, TIG-welded and heat-treated 6061-T6 aluminium tube set, which reduces the weight while at the same time providing great rigidity. Vital parts like the wheel forks, X-folding cross and handles are all made from high-grade aluminium.



Trionic Flex

A conventional 4-wheeled rollator always hits problems on uneven ground, as one of the four wheels sooner or later loses its contact with the surface.

The new "Trionic Flex" solution of the Walker 12er/14er has four wave springs integrated between the frame and the folding cross. This feature gives the Walker 12er & 14er an extra bit of torsional flex, so that it better remains standing on all four wheels.

Ergonomic benefits

When compared with other rollators the seat of the Walker is positioned further to the front, which means that you can walk closer to the grip handles and you automatically get a more upright and ergonomically correct body posture. The handles have height-indication marks and can be set at your exact preferred height, without any fixed increments, and they can also be turned to your chosen angle.

Transport

The Walker has detachable wheels with quick-releases. If required, you can minimize the transport weight and dimensions by quickly and easily detaching the wheels from the frame.

Wheel Sizes, Tires & Grips

Choose between multiple wheel sizes, tire types and grips in order to get the Walker custom-built according to your requirements and preferred activities. You can also choose

Walker Features

All Walker models include the following features:

- Trionic Syncro synchronized steering system.
- Screws, nuts, washers, axles and small parts made from stainless steel.
- Ackermann steering geometry.
- Detachable wheels with quick-releases.
- X-folding mechanism.
- Frame from machine-formed 6061-T6 aluminium tubing.
- Increment-free adjustment of handle height and angle.
- Adjustable grip handles (angle).
- Height adjustment through quick-release levers.
- All-terrain hub brakes for all weather conditions.
- All-aluminium brake levers with motion and parking brake function.
- Water-repellent and exchangeable textile accessories made from 600D polyester.
- 360° visible reflector for added safety when walking at night.

More product specific information with video clips is available at www.trionic.ca.



Trionic Walker

The Trionic Walker 9er, with 9" wheels (23 cm), is the most compact model and it is very well suited for both indoor- and outdoor use. Thanks to its hub brakes, the 9er is equipped with solid PU-tires of a much softer compound than other premium rollators. The tires are about 50% softer than the competition, and this of course raises your level of comfort when you for example hit obstacles or walk along cobblestone streets. They also have a distinct tread that offers you a secure grip and optimal safety.

With the Trionic Walker 12er you take a step up in wheel size, to 12" (31 cm), and this model is a great outdoor rollator. In order to raise your comfort level another notch, the Walker 12er is equipped with air-pressured tires. The 12er is your perfect companion for going grocery shopping, for taking daily strolls and for enjoying the outdoors. Depending on where you will mostly use your Walker 12er, you have the option to choose from 3 different tires.

By introducing the Trionic Walker 14er, Trionic takes yet another step into the future. It is the first ever rollator with 14" big wheels (36 cm). A bigger wheel simply performs better on uneven ground. With the Walker 14er you'll enjoy an unprecedented rollator experience.

On pages 24-25 you'll find several tire and grip options for your Walker.

For detailed technical specifications, measurements and weights, please see page 30.



Trionic Walker 9^{er} – the Combi Rollator

The Walker 9er is your perfect choice when you are looking for a rollator to use both indoors and outdoors. With its 9" wheels and extra-soft PU-tires you'll enjoy high comfort and a well-balanced performance. It's light and compact for easy storage and transport.

Size/Grip height:
User body height:
Wheel size:
Dimensions & weight:
Colours:
Including:





The new Walker 12er and Walker 14er are 2.6 lbs lighter



Trionic Walker 12er – the Outdoor Rollator

Are you looking for a rollator that copes with cobblestones and bumpy streets? Then the Walker 12er will be your faithful companion when you are heading outdoors. Its 12" wheels and air-filled tires guarantee a very comfortable ride, and it's perfect for every urban and outdoor area. Lightweight Trionic Primera tires with puncture protection. Stainless steel ball bearings, axles and small parts.

Sizes/Grip height:	S 28 - 33.9"	M/L	30.0-37.8"		
User body height:	S 4'8" - 5'8"	M/L	4'12"-6'4"		
Wheel size:	12" / 31 cm				
Dimensions & weight:	See page 30				
Colour size S:	Red with black/b	black/g	rey		
Colours size M/L:	Silver with black,	/black/	grey • navy/black/red		
Including:	Seat				
	Basket w. zip-pocket & 360° reflective band				



Trionic Walker 14^{er} – the Terminator

When you want nothing but the best. The Trionic Walker 14er is the world's first rollator with 14" big wheels. We call it the "terminator" as it leaves all other rollators trailing behind it. This machine makes you walk longer distances at a quicker pace and with greater ease! Lightweight Trionic Primera tires with puncture protection. Stainless steel ball bearings, axles and small parts.

Size/Grip height: User body height: Wheel size: Dimensions & weight: Colours: Including: M/L 32.0-39.4" M/L 5'4"-6'7" 14" / 36 cm See page 30 Silver with black/black/grey · green/black/red Seat Basket w. zip-pocket & 360° reflective band





What is a Veloped?

The Veloped is a very clever Swedish invention, best described as a "walking tricycle". It is tested and fully approved as a "rollator", but its design and its features differ a lot from conventional rollators. Through its unique design the Veloped promotes physical activity and outdoor activities. It is designed with a great focus on walking, and it is the first true off-road capable rollator/walker available in the world. When walking in uneven terrain, the Veloped beats any conventional rollator – hands down.

Off-Road 3-Wheel Design

The Veloped's 3-wheel design results in permanent ground contact on any uneven surface. Just like a 3-wheeled baby jogger with big wheels, the Veloped is a true all-terrain machine.

The Veloped is in fact the only walker/rollator on the market that is truly off-road capable.

All-Aluminium Design

The Veloped's framework is made from high-grade aluminium only, and it has no plastic parts whatsoever. The new frame is TIG-welded from hydroformed, butted and heat-treated 6061-T6 aluminium tubing, which keeps the weight low while at the same time offering great strength. Important parts like the fork, linkarm and wheel adjuster of the Trionic Climbing Wheel, the grip bar and the telescopic folding stay are all made from high-grade extruded and heat-treated aluminium.

Wheel Adjuster

The climbing wheel is adjustable for different conditions and terrain. For superior maneuverability, easy turning and with retained climbing capacity the front wheel is raised with a simple adjustment handle. In its other position, the all-terrain mode, the wheel has a suspension/climbing function with increased direction steadiness; perfect for the off-road walk.

Folding/Transport

Through a quick push of the button on the telescopic folding stay, you can fold your Veloped to pass narrow passages. To fully fold it for transport or storage, press the second button on the telescopic stay and fold the frame all the way. You can also push the grip bar to the bottom position and detach the wheels by using the quick release axles in order to minimize the transport dimensions.

Ergonomic Benefits

When the seat is slid forward the Veloped provides plenty of walking space within the frame. This means that you can walk closer to the grip handles, which automatically results in a more upright and ergonomically correct body posture. The U-shaped and ergonomic grip bar offers several grip positions, and while seated it acts as a back support. The grip bar can be freely height-adjusted in order to suit your body height perfectly.

General Veloped Features

All Veloped models include the following features:

- Trionic Climbing Wheel with climbing/suspension function.
- Dual-Mode adjustability for off-road or city use.
- 3-wheeler design for permanent ground contact.
- All-aluminium brake lever with motion and parking brake function.
- Detachable wheels with quick-releases.
- Frame from oversized 6061-T6 aluminium tubing.
- Ergonomic gripbar that doubles as a back rest.
- Height adjustment through quick-release levers.
- All-terrain hub brakes for all weather conditions.
- 2-step telescopic folding mechanism.
- Water-repellent and exchangeable textile accessories through the use of push buttons, velcro and/or buckles.
- Reflective front badge and reflective tire logos for added safety when walking at night.
- Screws, nuts, washers, axles and small parts made from stainless steel.

More product specific information with video clips is available at trionic.ca.







The Veloped Sport 12er offers you support and safety when you are exercising outdoors. It greatly reduces the impact and stress caused on your ankles, knees, hips and back and it'll make you last the distance. Equipped with 12" big wheels and air-pressured tires.

Sizes/Grip heights: User body height: Wheel size: Dimensions & weight: Colours: Including: M 29.5"-37.4" L 32.3"-40.2" M 4'1"-6'2" L 6'1"-6'9" 12" / 31 cm See page 31 navy/red/white · black/white/white 10" deep seat for maximum walking space · Basket w. zip-pocket & rain cover With its 14" big wheels the Veloped Sport 14er excels when you head for even rougher terrain, and it copes with higher obstacles than its little sibling. With this model you'll float over any uneven ground and it offers you unique possibilities to get or stay in shape.

Sizes/Grip heights: User body height: Wheel size: Dimensions & weight: Colours: Including: M 30.3"-38.2" L 33.1"-41.0" M 5'1"-6'4" L 6'1"-6'11" 14" / 36 cm See page 31 navy/red/white · black/white/white 25 cm deep seat for maximum walking space · Basket w. zip-pocket & rain cover

Veloped Sport

The Veloped Sport will be your faithful supporter when you practise the world's most underrated physical exercise – walking. Walking may be the closest thing to perfect exercise, and it's a great way to get in shape and to lose weight.

The intensity level of fitness walking is optimal for improving your health, while at the same time achieving a high fat burning effect. If your aim is to lose weight and improve your fitness at a steady rate, fitness walking is an excellent choice. You can exercise a little longer each time, and instead of wearing yourself down, like with more intensive workouts, you will enjoy faster progress when working towards your personal goal. For maximum health and fat burn effect, the ideal pulse rate is 55-70% of your maximum heart rate.

The Veloped Sport offers you comfort, support and safety while exercising. It greatly reduces the impact and stress on your body that can be caused by fitness walking. It has the same advantages as a mountain bike has to cyclists and you can head for green areas, walking paths, exercise trails and the forest. Now you can fully enjoy fitness walking – feeling safe and free to go wherever you want.

Vieloped warking

Get started today!

To get going, ask for our free 12-Week Walking Plan! You can also download the walking plan, get extensive exercise information and learnmore at www.trionic.ca.

Jette Dubois from Holland enjoying her daily exercise walk. FIFON



Veloped Tour

With the Veloped Tour you will enjoy daily comfort and practicality. It is very well suited for taking a stroll in the park, for going shopping and for sight-seeing. The roomy basket loads your groceries from the supermarket, and bringing them home is a joy with all the benefits of the climbing wheel and air-filled tires.

To provide optimal seating comfort, the seat depth has been increased to 12 inch. The longer seat, in combination with rigid carbon fiber inserts, offers plenty of support along the legs.

The colour design is sombre and elegant, with silver stripes and matching red lining. The stripes are made from 3M Scotchlite reflective fabric, and when walking at night the reflective stripes offer increased safety in traffic.

The Veloped Tour also comes with a practical rain cover that protects the basket and seat. The cover has a 2-colour design, and you can choose between silver and black depending on what side you turn outwards.



The Veloped Tour 12er is your perfect partner for everyday life. You can walk your dogs, go shopping and take a stroll in the local park. The 12" inch big wheels and air-filled tires give you comfort and security when walking along cobblestone streets, and it's easy to pass curbs and other obstacles along your path.

Sizes/Grip heights: User body height: Wheel size: Dimensions & weight: Colour: Including: M 29.5"-37.4" L 32.3"-40.2" M 4'1"-6'2" L 6'1"-6'9" 12" / 31 cm See page 31 black/reflex/red 12" deep comfort seat · Basket with zip-pocket and rain cover · Basket/seat rain cover By stepping up to 14" big wheels you'll enjoy even better driving characteristics. If you'll use your Veloped Tour both on- and off-road, the new Trionic X-Country tires are perfect. If you'll use it mainly on tarmac, please have a look at the optional Schwalbe Big Apple tires on page 25.

Sizes/Grip heights: User body height: Wheel size: Dimensions & weight: Colour: Including: M 30.3"-38.2" L 33.1"-41.0" M 51"-6'4" L 6'1"-6'11" 14"/36 cm See page 31 black/reflex/red 12" deep comfort seat · Basket with zip-pocket and rain cover · Basket/seat rain cover



The Veloped Trek 12er offers you comfort, support and safety during your outdoor adventures. It greatly reduces the impact and stress caused by trekking and hiking. With a Veloped Trek you'll enjoy the same benefits as when riding a mountain bike, and you can head for mountains, hiking trails, forests, dirt roads and green areas.

Sizes/Grip heights: User body height: Wheel size: Dimensions & weight: Colours: Including: M 29.5"-37.4" L 32.3"-40.2" M 4'1"-6'2" L 6'1"-6'9 12"/31 cm See page 31 green/black/orange · navy/black/yellow 10" deep seat · Basket with zip-pocket and rain cover · Seat bag · Backpack One of Trionic's biggest Veloped hiking fans is Hilary Ambrose from England. Her Veloped Trek goes under the nickname "Monster", and the new Veloped Trek 14er is indeed the Monster among walkers. It eats dirt roads and hiking trails for breakfast and will bring you anywhere you want to go.

Sizes/Grip heights: User body height: Wheel size: Dimensions & weight: Colours: Including: M 30.3"-38.2" L 33.1"-41.0" M 5'1"-6'4" L 6'1"-6'11" 14" / 36 cm See page 31 green/black/orange - navy/black/yellow 10" deep seat - Basket with zip-pocket and rain cover - Seat bag - Backpack

Veloped Trek

Load up your Veloped Trek and head for exciting adventures. It offers a multitude of options to carry all your gear when you are heading out to explore the world. With its versatile seat bag and backpack, the Veloped Trek is your perfect companion for hiking, camping, and travelling.

The seat bag snaps onto the seat, and is easily folded forward when you want to sit down for a break. It has three separate inner compartments, with plenty of room for your camera, maps, your valuables and a water bottle or thermos. The seat bag can also be carried as a shoulder bag, and has a padded detachable strap. One zippered outer compartment and two mesh pockets swallow your smaller items. A top elastic cord holds any extras.

The backpack fits nicely into the basket, and with a full-length zip you have easy access to the main compartment. The shoulder straps can be hidden in a separate pocket. Your MP3-player fits into the small top zippered pocket, and there is a sleeve for your head phones. In addition there is a zippered outer compartment and two mesh side pockets that fit two water bottles.

The seat cover, basket, seat bag and backpack are made from sturdy and water-repellent 600D polyester, with reflective lining for safety when walking in traffic at night. Alex Morley and his wife enjoying a hike in Oregon, USA.

Thomas having a rest while playing golf at Sigtunabygdens Golf Club, Sweden.



Veloped Golf

Stay in the game and do yourself a favour by walking the golf course instead of riding an electric golf cart. With the Veloped Golf you'll last the distance playing 18 holes while getting plenty of healthy exercise at the same time.

With all the benefits of the climbing wheel, the suspension, the 3-wheel design and its air-filled tires, the Veloped Golf finds itself very much at home on the golf course.

The Veloped Golf has two golf bags that are attached to the main frame. It loads a full set of 14 golf clubs, and the bags can either be closed or left open at the front depending on your preferred way of carrying your clubs. The basket loads your water bottle, some extra clothing and whatever you might need along the fairway.



Sizes/Grip heights:	12er M 12er L 14er M 14er L	29.5"-37.4" 32.3"-40.2" 30.3"-38.2" 33.1"-41.0"			
User body height:	12er M 12er L 14er M 14er L	4'1"-6'2" 6'1"-6'9" 5'1"-6'4" 6'1"-6'11"			
Wheel sizes:	12" 14"				
Dimensions & weight: Including: Colour:	See page 31 10" deep seat · Basket with zip-pocket & rain cover · Two detachable golf bags navy/black/red				



Sizes/Grip heights:	12er M 12er L 14er M 14er L	29.5"-37.4" 32.3"-40.2" 30.3"-38.2" 33.1"-41.0"			
User body height:	12er M 12er L 14er M 14er L	4'1"-6'2" 6'1"-6'9" 5'1"-6'4" 6'1"-6'11"			
Wheel sizes:	12" / 31 cm 14" / 36 cr				
Dimensions					
& weight:	See page 31 10" deep seat • Basket with zip-pocket & rain cover • Rifle bag • Backpack				
Including:					
Colour:	green/bla	ick/orange			

Veloped Jakt

The Veloped Jakt will be your trustworthy companion on your hunting trip. With its offroad capacity and climbing wheel, you can easily overcome roots, gravel, rock and dirt on your way to the hunting ground. The suspension, air-filled tires and the 3-wheel-design offer unbeatable comfort and consistent ground contact.

Even if you are a seasoned hunter, going on a hunting venture can be a major challenge if you suffer from knee problems, back pain, weakened joints, overweight or age-related issues.

Many hunting areas are not very easily accessed. You are required to hike along trails and dirt roads in order to get there. A proper "hunting walker" needs to perform well in all off-road terrain. Standard issue walkers and rollators are not at all intended for this kind of terrain, whereas the Veloped Jakt is. The rifle bag, that attaches to the main frame, relieves your back, and the seat is very practical when waiting for your prey.



You are unique

We are all different, and we have different preferences and requirements. The Velopeds & Walkers are all equipped with high-quality tires and very comfortable grips, and in addition Trionic offers you optional grip, tire and tube choices.

Grips

All Velopeds & Walkers are equipped with Ergon GP1 grip handles. They offer you a very comfortable grip while at the same time distributing the weight you put on them in an ergonomically correct way. The grip handles are available in two sizes (S for hand size 6.5-8.5 and L for hand size 8.5-12), and feature a forged aluminium clamp for quick and easy adjustment of your preferred grip angle.

Tires

The Velopeds are equipped with Trionic X-Country air-pressured tires as standard. This is a lightweight allround tire with a mountainbike-like tread that suits both off-road & on-road use. The X-Country has a puncture protection belt and reflective logos for added safety.

The Walker 12er & 14er come with the new Trionic Primera tire as standard. The Primera has a fine tread that suits tarmac, cobblestones and gravel, and thanks to its good flexibility it offers a very good comfort. It has a puncture protection belt and reflective logos.



Optimize your product



Grip Options

The grips of a walking aid are of major importance. Therefore we offer you optional grips of different sizes and made from different materials.

Ergon GP1 Grips

- German award-winning design.
- Available in two sizes, S for hand size 6.5-8.5 and L for hand size 8.5-12.
- Can be adjusted to your preferred angle.

Ergon GP1 BioKork Grips

The GP1 BioKork grips offer the same ultimate ergonomics for your hands as the standard GP1 grips. However, the unique qualities of cork mean that these grips are antibacterial and hypoallergenic.



Tyre Options

Depending on where you will mostly use your Veloped or Walker, you can choose between different tyres that suit different types of ground. All tyres offered are equipped with puncture protection.

Veloped Tyre Options

- Trionic X-Country 12"/14" Off-road/all-round
- Schwalbe Big Apple 12"/14" · City

Walker Tyre Options

- Trionic Primera 12"/14" · All-round/City
- Schwalbe Big Apple 12"/14" · City

Accessories & Parts

Multi-function Pump

This smart pump can be used both as a hand pump and a floor pump. A foldout foot pad holds the pump steady for better leverage. The T-handle provides a comfortable grip, and the flexible hose prevents the valve from being damaged or broken. Includes frame bracket.



Veloped Cane/Gear Holder

2-piece gear holder that also works as a safety strap when the Veloped is folded for transport. The gear holder is easily attached to the frame with Velcro straps.

Walker Cane Holder

2-piece cane holder that is attached to the semi-vertical frame tube of the Walker. The upper part has a smart quick-lock function. and the cane can quickly be attached and removed.

Veloped Storage

Rain Cover

Water-proof cover for storing your Veloped outdoors. Three hooks with elastic cords secure the cover so that it stays in place.



Colour: silver.

Veloped Seat/Basket

Rain Cover

A versatile rain cover that protects the seat and basket, with three options to attach the rear end.



2-colour design: silver/black.

Walker Seat/Basket

Rain Cover A detachable rain cover that protects the seat and basket of your Walker.



2-colour design: silver/black.



Veloped Bottle

Holder

Universal size plastic bottle holder that mounts to the frame or gripbar of the Veloped with a quick-release.

Colour: black.

Walker Bottle Holder

This aluminium bottle holder is open towards the side, so that it's easier to remove and enter the drinking bottle. The bottle holder can be adapted for rightor left-side attachment to the Walker frame.

Colour: black.

Sports Bottle

This sturdy water bottle has a screw-on cap with a lockable spout. The triangular shape of the cap guarantees a secure grip. The bottle is easy to fill and clean. Volume: 50 cl.

Colour: transparent/black/red



THEMONIE



LED Front Light

See and be seen with this bright LED light. The light offers two options, pulsing or steady light. Comes with a rubber band that offers multiple ways of attachment. Two versions available:

- Veloped LED front light w. 0° angle
- Walker LED front light w. 10° angle



LED Rear Lights

These two rear LED lights let you be seen everywhere. With flashing or steady mode, and with flexible rubber straps that enable multiple attachment positions. Two versions available:

- Veloped LED rear lights w. 0° angle
- Walker LED rear lights w. 10° angle



Both brakes (R+L) are operated with one single brake lever.The single-hand brake is intended for users who can operate the brakes with only one hand. For right- or left-side assembly.

At the order of a single-hand brake system, it is fully assembled to the Walker/Veloped before delivery.



Dimensions: Ø10 x L 900 mm



Accessories & Parts





Trionic Primera Tyre with Puncture Protection

Trionic X-Country Tyre

Dimension:

12 x 2.2" (310 x 57 mm)

14 x 2.2" (360 x 57 mm)

with Puncture Protection

Dimension: 12 x 2.00" (310 x 50 mm) 14 x 2.00" (360 x 50 mm)



Schwalbe Big Apple Tyre with Puncture Protection

Dimension: 12 x 2.00" (310 x 50 mm) 14 x 2.00" (360 x 50 mm)

Ergon GP1 Grips

Ergon's original and ground breaking GP1 grip continues to set standards. Offering optimal pressure distribution and available in two sizes, it features a forged aluminum clamp for a quick and easy adjustment of your preferred grip angle.

Grip size: S (Hand size: 6.5-8.5) L (Hand size: 8.5-12)

Ergon GP1 BioKork Grips

The GP1 BioKork offers the same ultimate ergonomics for the hand as the GP1 grip. However, the unique qualities of cork mean that the grip is antibacterial and hypoallergenic. Cork's natural damping qualities make it an ideal grip material.

Grip size:

S (Hand size: 6.5-8.5) L (Hand size: 8.5-12)



Veloped Trek

– Basket Bag/Backpack

20-litre backpack that fits nicely into the Trek basket. The padded shoulder straps and waist belt can be hidden in a rear side pocket. The bag has a large outer zip-pocket and



two side mesh pockets for water bottles. The top zip-pocket with sleeve fits your MP3-player and head phones.

Colours: green/black/orange navy/black/yellow

Veloped Trek – Seat Cover

and Basket

Seat cover and basket that match the colours of the backpack and seat bag. Seat cover with 4 push buttons for snapping on the seat bag. Seat depth: 10". The basket has a black rain cover integrated with the zip-pocket. Both made from water-repellent and durable 600D polyester.

Colours: green/black/orange navy/black/yellow



28



Veloped Trek – Seat/Shoulder Bag

10-litre bag that snaps onto the Veloped Trek seat with push buttons. It has a padded shoulder strap and can favourably be used as a shoulder bag. The main compartment can be divided into three separate parts to fit your camera, valuables and a water bottle or thermos. The outer zip-pocket and two mesh pockets allow quick access to smaller items.

Colours:

green/black/orange navy/black/yellow



Schwalbe Inner Tube

Sizes: 12", 14"

Veloped Tour – Comfort Seat Cover

and Basket Comfort seat with Carbon inserts for extra support. 3M Scotchlite reflective stripes. 30 cm seat depth. Roomy basket with a separate zip-pocket for your valuables. Both made from 600D polyester.



Basket colour: black/red

Veloped Sport – Seat Cover

and Basket

Seat cover that ensures full leg movement when fitness walking. Seat depth: 25 cm. Matching basket with zip-pocket. Perfect for some extra clothing when you are out exercising. Both made from durable 600D polyester with matching stripes.

Colours:

black/white/white navy/red/white





Veloped Golf – Seat Cover, Basket, Golf Bags & Backpack

Seat cover with 10" depth for maximum walking space. The basket has a rain cover integrated with the zip-pocket. The golf bag set consists of two pencil bags, each with two compartments.

Colour: navy/black/red

Walker Backrest

The ergonomic and padded backrest, made from water-repellent and durable 600D polyester, fits the rollators Walker 9er, 12er and 14er. It offers a comfortable and safe support when you are seated. The length/depth of the backrest is adjustable.





Trionic Walker - Technical Information

The Trionic Walker is tested and fully approved according to ISO 11199-2:2005, the relevant standard for rollators. It is registered as a "Class 1 Medical Device", and it bears a CE-mark. 10 Year Guarantee, parts subject to wear excluded.

Dimension Grip height	Walker 9er M 29.0-37.0"	Walker 12er S 28.0-33.9"	Walker 12er M/L 30.0-37.8"	Walker 14er M/L 32.0-39.4"	All-Terrain Brake	•	
User body height	4'10"-6'2"	4'8" - 5'8"	5'0"-6'4"	5'4"-6'7"			
v v		31.9"	31.9"	33.5"	X-Folding Mechanism	•	•
Length	26.8"						
Width	26.0"	27.2"	27.2"	27.2"	Trionic Syncro Steering	•	•
Grip width	18.5"	19.7"	19.7"	19.7"			
Seat height	24.4"	22.0"	24.4"	24.4"	Quick-Release Wheels	•	•
Seat width	16.9"	18.1"	18.1"	18.1"			
Wheel size	9"	12"	12"	14"	Seat	•	•
Dimensions folded					360° Reflector Band	•	•
Height	32.3"	33.9"	34.6"	35.4"	Basket with Zip-Pocket		
Length	26.8"	31.9"	31.9"	33.5"	Basket with Zip-Pocket	•	•
Width	11.0"	11.0"	11.0"	11.0"	Parking Brake	•	•
Transport dimensions -	folded without wheel	S			Motion Brake	•	•
Height	29.1"	29.9"	29.9"	29.9"			
Length	24.0"	26.4"	26.4"	26.4"	Air-Pressured Tires		•
Width	11.0"	11.0"	11.0"	11.0"	Extra-Soft PU-Tires	•	
Transport weight – with	outwheels						
		1/11/ba	1/ 7 lba	1/ Elles			
Transport weight	14.5 lbs	14.1 lbs	14.3 lbs	14.5 lbs			
Max user weight	330 lbs	330 lbs	330 lbs	330 lbs	NOTE: If you compare the weight	s of the	o Tric
Total weight	20.9 lbs	20.0 lbs	20.3 lbs	21.8 lbs	NOTE. If you compare the weight		

NOTE: If you compare the weights of the Trionic Walkers with other rollators, please note that the Trionic Walkers are equipped with high-performance hub brakes and premium tires.

Walker 14er Walker 12er Walker 9er

•

•

•

The tires and brakes of the Walkers offer you a very high comfort level and optimal driving characteristics; these features add some weight when compared with conventional, low cost rollator parts.

Trionic Veloped – Technical Information

The Trionic Veloped is tested and fully approved according to ISO 11199-2:2005, the relevant standard for rollators. It is registered as a "Class 1 Medical Device", and it bears a CE-mark. 10 Year Guarantee, parts subject to wear excluded.

Dimension	12er Medium	12er Large	14er Medium	14er Large	All-Terrain Brake System	•		•	•	•
Grip height	29.5-37.4" 4'11"-6'2"	32.3-40.2" 6'1"-6'9"	30.3-38.2" 5'1"-6'4"	33.1-41.0" 6'1"-6'11"		-				
User body height	411-62 40,1"	40.1"	42.1"	42.1"	2-step Folding Mechanism	•	•	•	•	•
Length Width	29.9"	31.1"	42.1 29.9"	42.1 31.1"						
Grip width	29.9 18.9"	18.9"	29.9 18.9"	18.9"	Climbing Wheel Technology	•	•	•	•	•
Seat height	22.8"	24.8"	23.6"	25.6"						
Seat width	17.5"	17.5"	17.5"	17.5"	Full Ground Contact Design	•	•	•	•	•
Wheel size	17.5	17.5	17.5	17.5	City/Off-road Mode	•			•	•
Wheel size	12	12	14	14			-	-		
Dimensions folded					Quick-Release Wheels	•	•	•	•	•
Height	29.5"	32.3"	30.3"	33.1"	Basic Seat, 25 cm depth	•		•	•	
Length	41.7"	41.7"	43.7"	43.7"	Basic Seat, 25 cm depth	—		<u> </u>		
Width	16.5"	16.5"	16.5"	16.5"	Comfort Seat, 30 cm depth		•			
Transport dimensions –	folded without wheel				Ergonomic Grip Bar	•	•	•	•	•
Height	23.2"	26.8"	23.2"	26.8"		-	-	-	-	
Length	33.1"	33.1"	33.1"	33.1"	Basket Rain Cover	•	•	•	•	•
Width	16.5"	16.5"	16.5"	16.5"	Seat/Basket Rain Cover		•			
						_				
Transport weight – with			167.0	17.6.1	Reflectors, 5 total	•	•	•	•	•
Transport weight	16.3 lbs	17.6 lbs	16.3 lbs	17.6 lbs	Reflective 3M Seat Stripes					
Max user weight	330 lbs	330 lbs	330 lbs	330 lbs						
Mux user weight	550 165	330 183	330 183	550 155	Basket with Zip-Pocket	•	•	•	•	•
Total weight*					Golf Bags			•		
Veloped Sport	25.1 lbs	26.4 lbs	26.2 lbs	27.6 lbs		_		-		
Veloped Tour	25.6 lbs	26.9 lbs	26.7 lbs	28.0 lbs	Rifle Bag					•
Veloped Trek	27.5 lbs	28.9 lbs	28.6 lbs	30.0 lbs		-				
Veloped Golf	30.8 lbs	32.2 lbs	32.0 lbs	33.3 lbs	Seat/Shoulder Bag				•	
Veloped Jakt	28.4 lbs	29.8 lbs	29.5 lbs	30.9 lbs	Basket Bag/Backpack				•	•
					Parking brake	•	•	•	•	•

Motion brake

Air-pressured tires

Veloped Colf Veloped Tour

• •

•

•

٠

.

•

•

•

Veloped Trek

Veloped Jakt

Veloped Sport



Trionic Canada c/o Dual FX Ltd.

692 Sumas Street Victoria, BC V8T 4S6 CANADA

Tel: 855-895-6635



www.trionic.ca distributor@trionic.ca

-

